

Netherlands Forensic Institute  
Ministry of Justice and Security

**PowerPlex Fusion 6C  
Profile analysis &  
interpretation**

*Corina Benschop*

12/11/2020  
Division Biological Traces

1

Netherlands Forensic Institute  
Ministry of Justice and Security

## Interpretation of mixed DNA profiles

- Clayton *et al.* (1998) Forensic Science International 91: 55–70

Step 1

Identify the presence of a mixture

↓

Step 2

Designate allele peaks

↓

Step 3

Identify the number of potential contributors

↓

Step 4

Estimate the relative ratio of the individuals contributing to the mixture

↓

Step 5


Consider all possible genotype combinations

↓


Step 6

Compare reference samples

2



3




4

2




## Unless...

- PCR inhibitors in DNA extract




May reduce PCR amplification efficiency for some or all loci

- Degraded DNA (low quality)




May make some allele targets unavailable (less intact target fragments)

- Low-template DNA (low quantity)



DNA-profile suffers from stochastic variation

- (Unbalanced) DNA mixtures and/or relatives involved



Alleles may be masked in the electropherograms

5


5



## PCR inhibitors

- PCR inhibitors interact with DNA or interfere with the DNA polymerase

100 pg/ $\mu$ l DNA007, 150 pmol/ $\mu$ l hematin:



*SGM+ 28c:*

*5  $\mu$ l input*

*2  $\mu$ l input*

*1  $\mu$ l input*

*0.5  $\mu$ l input*



6

6

3

## Prevent PCR inhibition

- Avoid unnecessary collection of inhibitors
- Decrease the PCR input
- Purify the DNA extract
- Use a different DNA polymerase
- Use an STR kit with optimised buffer  
e.g. NGM, PPF6C





7

7

## DNA degradation

- The DNA may be highly fragmented or chemically modified
- Different DNA degradation mechanisms
- Due to *e.g.* long time exposure to UV-light, high temperatures



## Possible solutions

- Multiple STR typing kits with different locus order
- mini-STRs, mtDNA, SNPs
- WGA prior to the PCR? Repair enzymes? Etc. ...


9

9



## Low template (LT) DNA

- LT DNA/ trace DNA/ touch DNA:  
*'any sample which falls below recommended thresholds at any stage of the process and can not be defined by a precise picogram amount'*
- DNA profiles with peaks below the stochastic threshold:  
*'threshold below which stochastic amplification artefacts can be expected'*




10

10



## On the stochastic threshold



\* Inception: VAL-00039 PowerPlex Fusion 6C sporen validatie

- Stochastic threshold of 800 rfus includes 98.9% of the single alleles

True genotype	Detected	Stochastic threshold risks
Homozygous {a,a}	1 allele {a}	<u>Too high</u> : allele is labeled {a,F} ⇒ risk of false inclusion
Heterozygous {a,b}	1 allele {a}	<u>Too low</u> : allele is labeled {a,a} ⇒ risk of false exclusion

11


11



## Heterozygote balance


Dependent on amount of DNA

High template



Balanced

Low-template



Peak imbalance  
can occur

Pg input	gemiddelde piekbalans per input	stdev	# loci
63pg (n=6)	75%	15%	110
125pg (n=6)	79%	15%	120
250pg (n=2)	85%	12%	40
500pg (n=2)	89%	9%	40
1000pg (n=2)	91%	7%	40

\* Inception: VAL-00039 PowerPlex Fusion 6C sporen validatie


12

12




## Stutter product formation


- Strand slippage during PCR



Normal replication






Increased -1 repeat  
unit (back stutter)



Increased +1 repeat  
unit (forward stutter)

Deletion caused by  
forward slippage

Insertion caused by  
backward slippage




13

13



## Replicate analysis




16

16



Netherlands Forensic Institute  
Ministry of Justice and Security


## Drop-outs in replicates



17

17






18

18



Next: profile interpretation



19

19




## Identify the presence of a mixture

- **By extra bands**
- **By peak imbalance**
- Exceptions: Genetic aberrations
  - Trisomy - three peaks of equal size
  - Somatic mutations - three peaks of unequal size
    - Can be tissue specific
    - Hairs have a higher rate
    - Many somatics will be missed as stutters
    - Not inherited (unless gametic)
  - Unbalanced heterozygote - e.g. XXY

*How often are genetic mutations observed?*  
**Trisomy** –  
[http://ed.mist.gov/biotech/rebase/tri\\_tab.htm](http://ed.mist.gov/biotech/rebase/tri_tab.htm)

Locus	Frequency
Amelogenin	0.1%
D21S11	0.0015%
D18S51	0.001%
DSS1179	0.004%
FGA	0.002%
VWA	0.0013%
TH01	0.00002%



20

20




## Identify the number of potential contributors


- Maximum allele count (MAC) approach:
  - At most 3 or 4 alleles/locus → at least two donors
  - At most 5 or 6 alleles/locus → at least three donors
  - >6 alleles/locus → the exact number of donors cannot be reliably determined
- Mischaracterisation rates with MAC can be high, specifically with high order mixtures ( $\geq 3$  contributors)
- The number of loci examined and their respective discriminating power, can reduce the mischaracterisation rates

21



21



22



23



G:\BI\BI\_DigiD

Excel file 'PPF6C\_TAC\_nC-tool\_v1'


24

24



## NOC-tool based on Machine Learning

- Features
  - Characteristics of the data
- Training, test and hold-out set
- Algorithms
- Model = algorithm + features
  - Performance



$$\sigma_{zz}) \quad (a-b)^2 = a^2 - 2ab + b^2 \quad \varepsilon_{zz} = \frac{\sigma_{zz}}{E} - \frac{\nu}{E} \sigma_{xx}$$
$$\frac{-4ac}{(x+a)^n} = \sum_{k=0}^n \binom{n}{k} x^k a^{n-k} \quad (a +$$

25

25



## Machine learning approach



Data collection



Features calculation



Feature selection



Parameter optimization



Model selection



Model validation


26

26




## Two NOC models

- PPF6C model (RFC19)
- Generic Model (2.0-Generic\_Model, RFC11)



27

27






28



## 2.0 Generic Model (kit independent)

### ESS core loci:

FGA, TH01, VWA, D1S1656, D2S441,  
D3S1358, D8S1179, D10S1248,  
D12S391, D18S51, D21S11,  
D22S1045.



29

29



## MAC vs PPF6C RFC19 model

True number of contributors (n per method)	MAC			PPF6C model		
	Under-estimated	Correctly estimated	Over-estimated	Under-estimated	Correctly estimated	Over-estimated
2 (n=90)	0%	67%	33%	0%	100%	0%
3 (n=88)	0%	97%	3%	3%	97%	0%
4 (n=89)	20%	77%	3%	17%	83%	0%
5 (n=87)	63%	37%	0%	40%	60%	not applicable
Total (n=354)	21%	69%	10%	15%	85%	0%

## MAC vs Generic RFC11 Model

True NOC	Percentage of correct predictions											
	PPF6C (NL)			NGMSE (AU)			GlobalFiler (FR)			GlobalFiler (SL)		
	n	MAC	Generic Model	n	MAC	Generic Model	n	MAC	Generic Model	n	MAC	Generic Model
Total	48	64.6%	95.8%	42	61.9%	71.4%	42	71.4%	78.6%	47	51.1%	70.2%
% increase correct predictions MAC vs Model	31.2%			9.5%			7.2%			19.1%		

31

31



## Implementation DNAXs

Number Of Contributors (NOC) estimate	≥5	<b>Run NOC Probability</b>	
Maximum Allele Count (MAC)	9		
Total Allele Count (TAC)	1	1A5.32 ≥5	0.658
Remaining alleles		Run 1 ≥5	0.658
Number Of Contributors (NOC) estimate	≥5	<b>Probabilities</b>	
Maximum Allele Count (MAC)	9		
Total Allele Count (TAC)	100		
Remaining alleles			
Max. remaining alleles per locus			
		1	0.000
		2	0.000
		3	0.010
		4	0.333
		≥5	0.658

Profile info

Derived from:

Complexity: Mixture

Duplicate of:

Profile quality: Complete

Sender: Unknown

DNA databank inclusion:

Number of Type 1 loci: 3

Number of Type 2 loci: 10

Number of Type 3 loci: 10

Number Of Contributors (NOC) estimate: ≥5

Maximum Allele Count (MAC): 9

Total Allele Count (TAC): 100

Remaining alleles:

Max. remaining alleles per locus:

Notes


32

32




**Step 3:**  
Estimate the relative  
ratio of the individuals  
contributing to the  
mixture, and, if  
possible deconvolve  
the mixture

## Estimate relative ratio of the contributors, if possible deconvolve



Indiv 1: 15/18

Indiv 2: 16/17



Indiv 1: 19/21

Indiv 2: 24/25

Or

Indiv 1: 19/24

Indiv 2: 21/25

Or

Indiv 1: 21/24

Indiv 2: 19/25

Or ...

33

33



## Deconvolution

Options for deconvolution/inference  
of a major contributor's profile, e.g.:

- Manual
- GeneMapper ID-X
- MasterMix
- **LoCIM-tool**
- Probabilistic using e.g.:
  - **EuroForMix**
  - DNAmixtures
  - STRmix
  - TrueAllele



Forensic Science International  
11 (1998) 41–53



Interpreting simple STR mixtures using allele peak  
areas

P. Gill<sup>a,\*</sup>, R. Sparkes<sup>a</sup>, R. Pinchin<sup>a</sup>, T. Clayton<sup>b</sup>, J. Whitaker<sup>b</sup>, J.  
Buckleton<sup>c</sup>

<sup>a</sup>Forensic Science Service, Priory House, Church Street North, Birmingham B190Q, UK  
<sup>b</sup>Forensic Science Service, Warburton Laboratory, Southack Way, Ashted Lane, Warburton, West Yorkshire  
LS23 4BW, UK

<sup>c</sup>ED Forensic, Mt Albert Science Centre, Private Bag 97-031, Auckland, New Zealand

Received 13 May 1997; received in revised form 9 October 1997; accepted 27 October 1997

Forensic Science International: Genetics 11 (2005) 154–165



Contents lists available at ScienceDirect  
Forensic Science International: Genetics  
journal homepage: [www.elsevier.com/locate/fgig](http://www.elsevier.com/locate/fgig)



LoCIM-tool: An expert's assistant for inferring the major contributor's  
alleles in mixed consensus DNA profiles

Corina C.G. Benschop<sup>a</sup>, Titia Sijen

<sup>a</sup>Department of Human Biological Sciences, Netherlands Forensic Institute, P.O. Box 24044, 2404 AK The Hague, The Netherlands




34

34




## LoCIM-tool

- **Locus Classification & Inference of the Major**
- Automated tool to deduce the major component's alleles in mixed (consensus) DNA profiles



35




36

Locus Classification

- Locus classification into three types:
    - Type 1
    - Type 2
    - Type 3
- (-Locus drop-out)

Increasing  
complexity  
for inferring  
the major



37

37

Locus Classification

- Locus classification based on:
  - PHs
  - Stochastic threshold
  - Heterozygote balance (HB)
  - Major to minor(s) ratio

> Theoretical PHs are assigned to the alleles in the consensus using the sum of the PHs in the individual amplifications.

	Alleles detected at locus D16S539				
	8	9	10	11	12
PHs amplification 1	123	484	584	66	62
PHs amplification 2	- <sup>1</sup>	690	705	75	-
PHs amplification 3	-	360	464	65	-
PHs amplification 4	135	725	582	65	-
Alleles in <i>n</i> /2 consensus	8	9	10	11	-
PHs assigned to alleles in the consensus (sum of PHs)	258	2259	2335	271	-

<sup>1</sup> The allele is not detected or below the detection threshold of 50 rfus.

38

38



Locus Classification - criteria

PPF6C 29c 1.2kV24s	PH alleles major	HB alleles major	Mhetero:m ratio	Mhomo:m ratio
Type 1	All ≥multiplied stochastic threshold	≥0.60	≥4:1	≥8:1
Type 2	-	≥0.60	≥2:1	≥4:1
Type 3	-	-	-	-

Increasing  
complexity  
for inferring  
the major  
↓

- Inferences include the allele with the highest summed PH (in the consensus) plus the alleles that have a PH within a specific percentage
  - Type 1 and type 2 loci: 50%
  - Type 3 loci: 33%


39

39



Typical examples

Type 1 loci



PPF6C 29c 1.2kV24s	PH alleles major	HB alleles major	Mhetero:m ratio	Mhomo:m ratio
Type 1	All ≥multiplied stochastic threshold	≥0.60	≥4:1	≥8:1
Type 2	-	≥0.60	≥2:1	≥4:1
Type 3	-	-	-	-


40

40



## Typical examples

### Type 2 loci



PPF6C 29c 1.2kV24s	PH alleles major	HB alleles major	Mhetero:m ratio	Mhomo:m ratio
Type 1	All $\geq$ multiplied stochastic threshold	$\geq 0.60$	$\geq 4:1$	$\geq 8:1$
Type 2	-	$\geq 0.60$	$\geq 2:1$	$\geq 4:1$
Type 3	-	-	-	-


41

41



## Typical examples

### Type 3 loci




PPF6C 29c 1.2kV24s	PH alleles major	HB alleles major	Mhetero:m ratio	Mhomo:m ratio
Type 1	All $\geq$ multiplied stochastic threshold	$\geq 0.60$	$\geq 4:1$	$\geq 8:1$
Type 2	-	$\geq 0.60$	$\geq 2:1$	$\geq 4:1$
Type 3	-	-	-	-

42

42

Inference of the Major




- Results for inference hold for consensus profiles generated from less amplifications and/or different CE setting, though with little shifts in the numbers for type 1, 2 and 3 loci
- LoCIM can be applied to single amplifications, though when an allele is not represented for the major the inference will be incorrect

43

43

Manual vs automated inference

- Results of the survey among 19 forensic scientists




44

44



## Locus classification and inference results related



45

45



## Boundaries of the LoCIM-tool?

- Test set of 131 DNA mixtures:
  - A) 2p HT mixtures in ratios 1:2 – 1:6
  - B) 2p LT mixtures in ratios 1:2 – 1:6
  - C) 2p mixtures (1:5, 1:10, 1:15) with relatives or degraded DNA
  - D) 3p and 4p degraded DNA mixtures having one major that is 2, 5 or 10 times more abundant than the other contributors
  - E) 2p, 3p and 4p mixtures having no major contributor
  - F) 3p and 4p mixtures representing multiple HT contributors
- Examined performance of the LoCIM-tool
- Explored whether the tool can handle more types of mixtures
- Examined consequences when applying the criteria to mixtures having no major and only LT contributors

46

46



## LoCIM test set - Inference results

- Type 1 loci: 99.9% correct inference  
0.1% missing allele (1/1003)
- Type 2 & 3 loci: Inference results as expected  
LT and degraded DNA: When alleles are not detected (in the consensus) they cannot be inferred

47

47



## Reliability of inferences

- Type 1 loci inferences are reliable
- Type 2 loci inferences are reliable if one can assume that:
  - 1) the mixture includes one major contributor
  - 2) the major contributor has no drop-out in the consensus
  - 3) the DNA profile shows no (severe) degradation
  - 4) PHs and balance are not (severely) affected by allele sharing
- Type 3 loci inferences may include extra alleles

**But also: Regard overall profile**

48

48



## Reliability of inferences → database actions?

- Type 1 loci inferences are reliable
  - *Can be stored*
- Type 2 loci inferences are reliable if one can assume that:
  - 1) the mixture includes one major contributor
  - 2) the major contributor has no drop-out in the consensus
  - 3) the DNA profile shows no (severe) degradation
  - 4) PHs and balance are not (severely) affected by allele sharing
  - *Can be stored if the above assumptions hold*
- Type 3 loci inferences may include extra alleles
  - *Advise against storage in the DDB, if based on LoCIM solely*


But also: Regard overall profile

49


49



## PPF6C LoCIM performance



Figuur 8. Resultaten van de 10 PPF6C profielen verkregen met de research tool. Het aantal loci per type is weergegeven alsook de correcte aflezingen (groen), hoofddonor allel(en) plus extra allel(en) afgeleid (geel), en incomplete aflezingen (oranje).




Figuur 9. Overzicht van de resultaten van de 10 PPF6C mengsels verkregen met de research tool. 1 = 20 pg. De bovenste grafiek geeft het soort locus typeringen per mengsel aan (type 1, type 2 en type 3) en de onderste grafiek het aantal correcte aflezingen (groen), hoofddonor allel(en) plus extra allel(en) afgeleid (geel), en incomplete aflezingen (oranje) per mengsel.


50

50





52



53

25



## Next: Read and practice

To read in advance:

- QOL-00434 Validatierapport PPF6C sporen
- QOL-00711 Interpretatie van autosomale STR DNA-profielen (t/m 7. Deconvolutie functionaliteit in DNAXs)

To read after this training:

- Paper PPF6C LT
- Paper LoCIM-tool
- Paper en validatierapport NOC-tool

And have fun!

